Nombre del documento: P21-094 Release Notes Hybrid Recon 5.0.0 Rev.1_SP

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Referencia FSN: P21-092-02 Field Safety Notice FSN0066683

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Hybrid Recon 5.0.0

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The screenshot shows the 'Program parameters' dialog box with the 'Recon' tab selected. The dialog is organized into several sections for different medical protocols:

- Cardiology protocols:** Includes fields for 'Study 1', 'Study 2', and 'Study 3'. It has input fields for 'Primary non-gated', 'Primary gated', and 'Secondary non-gated'. A radio button for 'Secondary NG' is selected. A 'Count limit [%]' spinner is set to 20. There are radio buttons for 'Outline' and 'Extend transverse FOV'.
- Lung protocols:** Includes input fields for 'Study 1' and 'Study 2'.
- Neurology protocols:** Includes input fields for 'Primary', 'Secondary', and 'Tertiary' (set to 'hermes_rec_default'). Radio buttons for 'Perform secondary reconstruction' and 'Perform tertiary reconstruction' are present.
- Oncology protocols:** Includes input fields for 'Primary' (set to 'onco_rec_default'), 'Secondary' (set to 'onco_nac_rec_default'), and 'Tertiary' (set to 'hermes_rec_default'). Radio buttons for 'Perform secondary reconstruction', 'Perform tertiary reconstruction', and 'Reconstruct whole FOV in NAC' are present.
- Mumap display zoom [%]:** Set to 100. **Uniform mumap display zoom [%]:** Set to 265.
- Warning options:**
 - Skip subset warning
 - Skip frame of reference warning (highlighted with a red box)
 - Reduce number of reconstruction methods
 - Skip energy window warning
 - Skip mumap truncation warning
- Mumap transformation:** Radio buttons for 'Automatic full 6 parameter', 'Automatic translation only' (selected), and 'Manual'.
- Mumap page screen cap label:** An empty text input field.

A 'Save' button is located at the bottom center of the dialog.

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Hermes Medical Solutions Ltd
Hermes Medical Solutions, Inc
Hermes Medical Solutions Canada, Inc
Hermes Medical Solutions Germany GmbH