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The screenshot shows the 'Program parameters' dialog box with the 'Recon' tab selected. The dialog is organized into several sections for different medical protocols:

- Cardiology protocols:** Includes fields for 'Study 1', 'Study 2', and 'Study 3'. It has input fields for 'Primary non-gated', 'Primary gated', and 'Secondary non-gated'. A radio button for 'Secondary NG' is selected. A 'Count limit [%]' is set to 20. There are radio buttons for 'Outline' and 'Extend transverse FOV'.
- Lung protocols:** Includes input fields for 'Study 1' and 'Study 2'.
- Neurology protocols:** Includes input fields for 'Primary', 'Secondary', and 'Tertiary' (set to 'hermes_rec_default'). Radio buttons for 'Perform secondary reconstruction' and 'Perform tertiary reconstruction' are present.
- Oncology protocols:** Includes input fields for 'Primary' (set to 'onco_rec_default'), 'Secondary' (set to 'onco_nac_rec_default'), and 'Tertiary' (set to 'hermes_rec_default'). Radio buttons for 'Perform secondary reconstruction', 'Perform tertiary reconstruction', and 'Reconstruct whole FOV in NAC' are present.
- Mumap display zoom [%]:** Set to 100. **Uniform mumap display zoom [%]:** Set to 265.
- Warning options:** A list of checkboxes: 'Skip subset warning', 'Skip energy window warning', 'Skip frame of reference warning' (highlighted with a red box), and 'Skip mumap truncation warning'. There is also an option to 'Reduce number of reconstruction methods'.
- Mumap transformation:** Radio buttons for 'Automatic full 6 parameter', 'Automatic translation only' (selected), and 'Manual'.
- Mumap page screen cap label:** An empty input field.

A 'Save' button is located at the bottom center of the dialog.

Indstillingen af alternativknappen "Skip Frame of Reference warning" [spring referencerammeadvarsel over] skal kontrolleres for at sikre, at den er indstillet til "off" [fra] for hver enkelt arbejdsgang – onkologi, neurologi, kardiologi og lunge.

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USA

Tlf.: +1 (866) 437-6372

Hermes Medical Solutions Canada, Inc

1155, René-Lévesque O., Suite 2500

Montréal (QC) H3B 2K4

Canada

Tlf.: +1 (877) 666-5675

Fax: +1 (514) 288-1430

Hermes Medical Solutions Germany GmbH

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48282 Emsdetten

Tyskland

Tlf.: +46 (0)819 03 25

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Hermes Medical Solutions Ltd
Hermes Medical Solutions, Inc